SCA-TEST®

versione 2.0

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The screenshot shows the SCA-TEST software interface. At the top, there is a blue header with the text "SCA-TEST®" and "Integrato 1° e 2° trimestre". Below the header, there is a navigation bar with a back arrow, four status indicators (one red, three grey), and a printer icon labeled "Integrato 1° e 2°".

Below the navigation bar, there is a section labeled "T21" containing a table of data and a 3D bar chart.

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The 3D bar chart shows two bars: a red bar for "ETA'" and a black bar for "INT". The y-axis ranges from 0 to 0.5. A horizontal green line is drawn at the 0.4 level.

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Bibliografia

1. Wright D, Spencer K, Kagan K K, Tørring N, Petersen OB, Christou A, Kallikas J, Nicolaides KH. First-trimester combined screening for trisomy 21 at 7-14 weeks' gestation. *Ultrasound Obstet Gynecol.* 2010 Oct;36(4):404-11.
2. Wald NJ, Rodeck C, Hackshaw AK, Walters J, Chitty L, Mackinson AM. First and second trimester antenatal screening for Down's syndrome: the results of the Serum, Urine and Ultrasound Screening Study (SURUSS). *J Med Screen.* 2003;10(2):56-104. Erratum in: *J Med Screen.* 2006;13(1):51-2.
3. Wapner R, Thom E, Simpson JL, Pergament E, Silver R, Filkins K, Platt L, Mahoney M, Johnson A, Hogge WA, Wilson RD, Mohide P, Hershey D, Krantz D, Zachary J, Snijders R, Greene N, Sabbagha R, MacGregor S, Hill L, Gagnon A, Hallahan T, Jackson L; First Trimester Maternal Serum Biochemistry and Fetal Nuchal Translucency Screening (BUN) Study Group. First-trimester screening for trisomies 21 and 18. *N Engl J Med.* 2003 Oct 9;349(15):1405-13.
4. Wald NJ, Rodeck C, Hackshaw AK, Walters J, Chitty L, Mackinson AM; SURUSS Research Group. First and second trimester antenatal screening for Down's syndrome: the results of the Serum, Urine and Ultrasound Screening Study (SURUSS). *Health Technol Assess.* 2003;7(11):1-77
5. Bindra R, Heath V, Liao A, Spencer K, Nicolaides KH. One-stop clinic for assessment of risk for trisomy 21 at 11-14 weeks: a prospective study of 15 030 pregnancies. *Ultrasound Obstet Gynecol.* 2002 Sep;20(3):219-25.
6. Schuchter K, Hafner E, Stangl G, Metzenbauer M, Höfingler D, Philipp K. The first trimester 'combined test' for the detection of Down syndrome pregnancies in 4939 unselected pregnancies. *Prenat Diagn.* 2002 Mar;22(3):211-5.
7. Niemimaa M, Suonpää M, Perheentupa A, Seppälä M, Heinonen S, Laitinen P, Ruokonen A, Ryyänänen M. Evaluation of first trimester maternal serum and ultrasound screening for Down's syndrome in Eastern and Northern Finland. *Eur J Hum Genet.* 2001 Jun;9(6):404-8.
8. Krantz DA, Hallahan TW, Orlandi F, Buchanan P, Larsen JW Jr, Macri JN. First-trimester Down syndrome screening using dried blood biochemistry and nuchal translucency. *Obstet Gynecol.* 2000 Aug;96(2):207-13.
9. Cicero S, Bindra R, Rembouskos G, Spencer K, Nicolaides KH. Integrated ultrasound and biochemical screening for trisomy 21 using fetal nuchal translucency, absent fetal nasal bone, free beta-hCG and PAPP-A at 11 to 14 weeks. *Prenat Diagn.* 2003 Apr;23(4):306-10.
10. Maiz N, Valencia C, Kagan KO, Wright D, Nicolaides KH. Ductus venosus Doppler in screening for trisomies 21, 18 and 13 and Turner syndrome at 11-13 weeks of gestation. *Ultrasound Obstet Gynecol.* 2009 May;33(5):512-7.
11. Kagan KO, Valencia C, Livanos P, Wright D, Nicolaides KH. Tricuspid regurgitation in screening for trisomies 21, 18 and 13 and Turner syndrome at 11+0 to 13+6 weeks of gestation. *Ultrasound Obstet Gynecol.* 2009 Jan;33(1):18-22.
12. Sonek J, Borenstein M, Dagklis T, Persico N, Nicolaides KH. Frontomaxillary facial angle in fetuses with trisomy 21 at 11-13(6) weeks. *Am J Obstet Gynecol.* 2007 Mar;196(3):271.e1-4.
13. Chaoui R, Benoit B, Mitkowska-Wozniak H, Heling KS, Nicolaides KH. Assessment of intracranial translucency (IT) in the detection of spina bifida at the 11-13-week scan. *Ultrasound Obstet Gynecol* 2009; 34: 249-252.
14. Rembouskos G, Passamonti U, De Robertis V, Tempesta A, Campobasso G, Volpe G, Gentile M, Volpe P. Aberrant right subclavian artery (ARSA) in unselected population at first and second trimester ultrasonography. *Prenat Diagn.* 2012 Jul 31:1-8.
15. Ozkaya O, Sezik M, Ozbasar D, Kaya H. Abnormal ductus venosus flow and tricuspid regurgitation at 11-14 weeks' gestation have high positive predictive values for increased risk in first-trimester combined screening test: results of a pilot study. *Taiwan J Obstet Gynecol.* 2010 Jun;49(2):145-50.
16. Bethune M. Literature review and suggested protocol for managing ultrasound soft markers for Down syndrome: thickened nuchal fold, echogenic bowel, shortened femur, shortened humerus, pyelectasis and absent or hypoplastic nasal bone. *Australas Radiol.* 2007 Jun;51(3):218-25
17. Nicolaides KH. Screening for fetal aneuploidies at 11 to 13 weeks. *Prenat Diagn.* 2011 Jan;31(1):7-15. doi: 10.1002/pd.2637.
18. Valinen Y, Marttala J, Sahravand M, Jarvela I, Honkasalo T, Laitinen P, Ryyänänen M. Adding ADAM12 in risk calculation program does not improve the detection rate of trisomies 18 and 13 in first trimester screening. *J Matern Fetal Neonatal Med.* 2012 Feb 6

19. Valinen Y, Laitinen P, Ranta J, Ignatius J, Jarvela I, Ryyänänen M. Effect of a new marker, ADAM12, on Down risk figures in first trimester screening. *J Matern Fetal Neonatal Med.* 2009 Jul;22(7):602-7.
20. Kublickas M, Crossley J, Aitken D. Screening for Down's syndrome in the first trimester: combined risk calculation, methodology, and validation of a web-based system. *Acta Obstet Gynecol Scand.* 2009;88(6):635-8.
21. H.P. Robinson. Sonar measurement of fetal crown-rump length as means of assessing maturity in first trimester of pregnancy. *Brit. Med. J.*, 4:28-31, 1973.
22. Salomon LJ, Chalouhi GE, Bernard JP, Ville Y; Société française pour l'amélioration des pratiques échographiques (SFAPE). Nuchal translucency thickness at 11-14 weeks of gestation: French charts and equations. *J Gynecol Obstet Biol Reprod (Paris).* 2009 Dec;38(8):635-41.
23. Spencer K, Nicolaides KH. Screening for trisomy 21 in twins using first trimester ultrasound and maternal serum biochemistry in a one-stop clinic: a review of three years experience. *BJOG.* 2003;110:27680.
24. Spencer K, Bindra R, Cacho AM, Nicolaides KH. The impact of correcting for smoking status when screening for chromosomal anomalies using maternal serum biochemistry and fetal nuchal translucency thickness in the first trimester of pregnancy. *Prenat Diagn.* 2004;24:169-73.
25. Liao AW, Heath V, Kametas N, Spencer K, Nicolaides KH. First-trimester screening for Trisomy 21 in singleton pregnancies achieved by assisted reproduction. *Hum Reprod.* 2001;16:1501-4.
26. Cuckle HS, Wald NJ, Thompson SG. Estimating a woman's risk of having a pregnancy associated with Down's syndrome using her age and serum alpha-fetoprotein level. *Br J Obstet Gynaecol.* 1987;94:38740.
27. Snijders RJM, Sundberg K, Holzgreve W, Henry G, Nicolaides KH. Maternal age- and gestation-specific risk for Trisomy 21. *Ultrasound Obstet Gynecol.* 1999;13:167-70.
28. Nicolaides KH. Patient-specific risk for chromosomal defects. In: Nicolaides KH, editor. *The 11136 weeks scan.* London: Fetal Medicine Foundation, 2004. pp. 1321.
29. Wald NJ. Biochemical screening. In: Brock DJH, Rodeck CH, Ferguson-Smith MA, editors. *Prenatal diagnosis and screening.* Oxford: Churchill Livingstone, 1992. pp. 563-78.
30. Nicolaides KH. Patient-specific risk for chromosomal defects. In: Nicolaides KH, editor. *The 11136 weeks scan.* London: Fetal Medicine Foundation, 2004. pp. 13-21.
31. Aitken DA, Crossley JA, Spencer K. Prenatal screening for neural tube defects and aneuploidy. In: Rimoin DL, Connor JM, Pyeritz RE, Korf BR, editors. *Principles and practice of medical genetics.* Churchill Livingstone, 2007. pp. 636-78.
32. Spencer K, Ong CY, Liao AW, Nicolaides KH. The influence of ethnic origin on first trimester biochemical markers of chromosomal abnormalities. *Prenat Diagn.* 2000 Jun;20(6):491-4.
33. Sahota DS, Leung TY, Fung TY, Chan LW, Law LW, Lau TK. Medians and correction factors for biochemical and ultrasound markers in Chinese women undergoing first-trimester screening for trisomy 21. *Ultrasound Obstet Gynecol.* 2009 Apr;33(4):387-93.
34. Spencer K, Cowans NJ, Spencer CE, Achillea N. A re-evaluation of the influence of maternal insulin-dependent diabetes on fetal nuchal translucency thickness and first-trimester maternal serum biochemical markers of aneuploidy. *Prenat Diagn.* 2010 Oct;30(10):937-40.
35. Faiola S, Tsoi E, Huggon IC, Allan LD, Nicolaides KH. Likelihood ratio for trisomy 21 in fetuses with tricuspid regurgitation at the 11 to 13 + 6-week scan. *Ultrasound Obstet Gynecol.* 2005 Jul;26(1):22-7.
36. Cicero S, Rembouskos G, Vandecruys H, Hogg M, Nicolaides KH. Likelihood ratio for trisomy 21 in fetuses with absent nasal bone at the 11-14-week scan. *Ultrasound Obstet Gynecol.* 2004 Mar;23(3):218-23.
37. Vranken G, Reynolds T, Van Nueten J. Medians for second-trimester maternal serum markers: geographical differences and variation caused by median multiples-of-median equations. *J Clin Pathol.* 2006 Jun;59(6):639-44.
38. Wald NJ, Densem JW, Smith D, Klee GG. Four-marker serum screening for Down's syndrome. *Prenat Diagn.* 1994 Aug;14(8):707-16.
39. Reynolds T. The triple test as a screening technique for Down syndrome: reliability and relevance. *Int J Womens Health.* 2010 Aug 9;2:83-8.

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Bibliografia

1. Goldenberg RL, Davis RO, Cutter GR, Hoffman HJ, Brumfield CG, Foster JM. Prematurity, postdates, and growth retardation: the influence of use of ultrasonography on reported gestational age. *Am J Obstet Gynecol.* 1989 Feb;160(2):462-70.
2. Gardosi J, Vanner T, Francis A. Gestational age and induction of labour for prolonged pregnancy. *Br J Obstet Gynaecol.* 1997 Jul;104(7):792-7.
3. Salomon LJ, Bernard M, Amarsy R, Bernard JP, Ville Y. The impact of crown-rump length measurement error on combined Down syndrome screening: a simulation study. *Ultrasound Obstet Gynecol.* 2009 May;33(5):506-11.
4. Gardosi J, Mongelli M. Risk assessment adjusted for gestational age in maternal serum screening for Down's syndrome. *BMJ.* 1993 Jun 5;306(6891):1509-11.
5. Verburg BO, Steegers EA, De Ridder M, Snijders RJ, Smith E, Hofman A, Moll HA, Jaddoe VW, Witteman JC. New charts for ultrasound dating of pregnancy and assessment of fetal growth: longitudinal data from a population-based cohort study. *Ultrasound Obstet Gynecol.* 2008 Apr;31(4):388-96.
6. Blaas HG, Eik-Nes SH, Bremnes JB. The growth of the human embryo. A longitudinal biometric assessment from 7 to 12 weeks of gestation. *Ultrasound Obstet Gynecol.* 1998 Nov;12(5):346-54.
7. Morin I, Morin L, Zhang X, Platt RW, Blondel B, Bréart G, Usher R, Kramer MS. Determinants and consequences of discrepancies in menstrual and ultrasonographic gestational age estimates. *BJOG.* 2005 Feb;112(2):145-52.
8. Wilcox M, Gardosi J, Mongelli M, Ray C, Johnson I. Birth weight from pregnancies dated by ultrasonography in a multicultural British population. *BMJ.* 1993 Sep 4;307(6904):588-91.
9. Papaioannou GI, Syngelaki A, Poon LC, Ross JA, Nicolaides KH. Normal ranges of embryonic length, embryonic heart rate, gestational sac diameter and yolk sac diameter at 6-10 weeks. *Fetal Diagn Ther.* 2010;28(4):207-19.
10. Pexsters A, Daemen A, Bottomley C, Van Schoubroeck D, De Catte L, De Moor B, D'Hooghe T, Lees C, Timmerman D, Bourne T. New crown-rump length curve based on over 3500 pregnancies. *Ultrasound Obstet Gynecol.* 2010 Jun;35(6):650-5.
11. Sladkevicius P, Saltvedt S, Almström H, Kublickas M, Grunewald C, Valentin L. Ultrasound dating at 12-14 weeks of gestation. A prospective cross-validation of established dating formulae in in-vitro fertilized pregnancies. *Ultrasound Obstet Gynecol.* 2005 Oct;26(5):504-11.
12. Wisser J, Dirschedl P, Krone S. Estimation of gestational age by transvaginal sonographic measurement of greatest embryonic length in dated human embryos. *Ultrasound Obstet Gynecol.* 1994 Nov 1;4(6):457-62.
13. Grisolia G, Milano K, Pilu G, Banzi C, David C, Gabrielli S, Rizzo N, Morandi R, Bovicelli L. Biometry of early pregnancy with transvaginal sonography. *Ultrasound Obstet Gynecol.* 1993 Nov 1;3(6):403-11.
14. Hadlock FP, Shah YP, Kanon DJ, Lindsey JV. Fetal crown-rump length: reevaluation of relation to menstrual age (5-18 weeks) with high-resolution real-time US. *Radiology.* 1992 Feb;182(2):501-5.
15. Goldstein I, Zimmer EA, Tamir A, Peretz BA, Paldi E. Evaluation of normal gestational sac growth: appearance of embryonic heartbeat and embryo body movements using the transvaginal technique. *Obstet Gynecol.* 1991 Jun;77(6):885-8.
16. Izquierdo LA, Kushnir O, Smith JF, Gilson GJ, Chatterjee MS, Qualls C, Curet LB. Evaluation of fetal sonographic measurements in the first trimester by transvaginal sonography. *Gynecol Obstet Invest.* 1991;32(4):206-9.
17. Silva PD, Mahairas G, Schaper AM, Schauburger CW. Early crown-rump length. A good predictor of gestational age. *J Reprod Med.* 1990 Jun;35(6):641-4.
18. Koornstra G, Wattel E, Exalto N. Crown-rump length measurements revisited. *Eur J Obstet Gynecol Reprod Biol.* 1990 May-Jun;35(2-3):131-8.
19. Vollebergh JH, Jongsma HW, van Dongen PW. The accuracy of ultrasonic measurement of fetal crown-rump length. *Eur J Obstet Gynecol Reprod Biol.* 1989 Mar;30(3):253-6.
20. Rossavik IK, Torjusen GO, Gibbons WE. Conceptual age and ultrasound measurements of gestational sac and crown-rump length in in vitro fertilization pregnancies. *Fertil Steril.* 1988 Jun;49(6):1012-7.
21. MacGregor SN, Tamura RK, Sabbagha RE, Minogue JP, Gibson ME, Hoffman DI. Underestimation of gestational age by conventional crown-rump length dating curves. *Obstet Gynecol.* 1987 Sep;70(3 Pt 1):344-8.

22. Selbing A, Fjällbrant B. Accuracy of conceptual age estimation from fetal crown-rump length. *J Clin Ultrasound*. 1984 Jul-Aug;12(6):343-6.
23. Selbing A. Gestational age and ultrasonic measurement of gestational sac, crown-rump length and biparietal diameter during first 15 weeks of pregnancy. *Acta Obstet Gynecol Scand*. 1982;61(3):233-5.
24. Nelson LH. Comparison of methods for determining crown-rump measurement by real-time ultrasound. *J Clin Ultrasound*. 1981 Feb;9(2):67-70.
25. Drumm JE, Clinch J, MacKenzie G. The ultrasonic measurement of fetal crown-rump length as a method of assessing gestational age. *Br J Obstet Gynaecol*. 1976 Jun;83(6):417-21.
26. Robinson HP, Fleming JE. A critical evaluation of sonar "crown-rump length" measurements. *Br J Obstet Gynaecol*. 1975 Sep;82(9):702-10.
27. Waldenström U, Axelsson O, Nilsson S. A comparison of the ability of a sonographically measured biparietal diameter and the last menstrual period to predict the spontaneous onset of labor. *Obstet Gynecol*. 1990 Sep;76(3 Pt 1):336-8.
28. Chervenak FA, Skupski DW, Romero R, Myers MK, Smith-Levitin M, Rosenwaks Z, Thaler HT. How accurate is fetal biometry in the assessment of fetal age? *Am J Obstet Gynecol*. 1998 Apr;178(4):678-87.
29. Taipale P, Hiilesmaa V. Predicting delivery date by ultrasound and last menstrual period in early gestation. *Obstet Gynecol*. 2001 Feb;97(2):189-94.
30. Chervenak FA, Skupski DW, Romero R, Myers MK, Smith-Levitin M, Rosenwaks Z, Thaler HT. How accurate is fetal biometry in the assessment of fetal age? *Am J Obstet Gynecol*. 1998 Apr;178(4):678-87.
31. Altman DG, Chitty LS. New charts for ultrasound dating of pregnancy. *Ultrasound Obstet Gynecol*. 1997 Sep;10(3):174-91.
32. Mul T, Mongelli M, Gardosi J. A comparative analysis of second-trimester ultrasound dating formulae in pregnancies conceived with artificial reproductive techniques. *Ultrasound Obstet Gynecol*. 1996 Dec;8(6):397-402.
33. Persson PH, Weldner BM. Reliability of ultrasound fetometry in estimating gestational age in the second trimester. *Acta Obstet Gynecol Scand*. 1986;65(5):481-3.
34. Selbing A, Kjessler B. Conceptual dating by ultrasonic measurement of the fetal biparietal diameter in early pregnancy. *Acta Obstet Gynecol Scand*. 1985;64(7):593-7.
35. Hadlock FP, Deter RL, Harrist RB, Park SK. Fetal biparietal diameter: rational choice of plane of section for sonographic measurement. *AJR Am J Roentgenol*. 1982 May;138(5):871-4.
36. Hadlock FP, Deter RL, Harrist RB, Park SK. Estimating fetal age: computer-assisted analysis of multiple fetal growth parameters. *Radiology*. 1984 Aug;152(2):497-501.
37. Chervenak FA, Skupski DW, Romero R, Myers MK, Smith-Levitin M, Rosenwaks Z, Thaler HT. How accurate is fetal biometry in the assessment of fetal age? *Am J Obstet Gynecol*. 1998 Apr;178(4):678-87.

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Bibliografia

1. Wright D, Spencer K, Kagan K K, Tørring N, Petersen OB, Christou A, Kallikas J, Nicolaides KH. First-trimester combined screening for trisomy 21 at 7-14 weeks' gestation. *Ultrasound Obstet Gynecol.* 2010 Oct;36(4):404-11.
2. Wald NJ, Rodeck C, Hackshaw AK, Walters J, Chitty L, Mackinson AM. First and second trimester antenatal screening for Down's syndrome: the results of the Serum, Urine and Ultrasound Screening Study (SURUSS). *J Med Screen.* 2003;10(2):56-104. Erratum in: *J Med Screen.* 2006;13(1):51-2.
3. Wapner R, Thom E, Simpson JL, Pergament E, Silver R, Filkins K, Platt L, Mahoney M, Johnson A, Hogge WA, Wilson RD, Mohide P, Hershey D, Krantz D, Zachary J, Snijders R, Greene N, Sabbagha R, MacGregor S, Hill L, Gagnon A, Hallahan T, Jackson L; First Trimester Maternal Serum Biochemistry and Fetal Nuchal Translucency Screening (BUN) Study Group. First-trimester screening for trisomies 21 and 18. *N Engl J Med.* 2003 Oct 9;349(15):1405-13.
4. Wald NJ, Rodeck C, Hackshaw AK, Walters J, Chitty L, Mackinson AM; SURUSS Research Group. First and second trimester antenatal screening for Down's syndrome: the results of the Serum, Urine and Ultrasound Screening Study (SURUSS). *Health Technol Assess.* 2003;7(11):1-77
5. Bindra R, Heath V, Liao A, Spencer K, Nicolaides KH. One-stop clinic for assessment of risk for trisomy 21 at 11-14 weeks: a prospective study of 15 030 pregnancies. *Ultrasound Obstet Gynecol.* 2002 Sep;20(3):219-25.
6. Schuchter K, Hafner E, Stangl G, Metzenbauer M, Höfingler D, Philipp K. The first trimester 'combined test' for the detection of Down syndrome pregnancies in 4939 unselected pregnancies. *Prenat Diagn.* 2002 Mar;22(3):211-5.
7. Niemimaa M, Suonpää M, Perheentupa A, Seppälä M, Heinonen S, Laitinen P, Ruokonen A, Ryyänen M. Evaluation of first trimester maternal serum and ultrasound screening for Down's syndrome in Eastern and Northern Finland. *Eur J Hum Genet.* 2001 Jun;9(6):404-8.
8. Krantz DA, Hallahan TW, Orlandi F, Buchanan P, Larsen JW Jr, Macri JN. First-trimester Down syndrome screening using dried blood biochemistry and nuchal translucency. *Obstet Gynecol.* 2000 Aug;96(2):207-13.
9. Cicero S, Bindra R, Rembouskos G, Spencer K, Nicolaides KH. Integrated ultrasound and biochemical screening for trisomy 21 using fetal nuchal translucency, absent fetal nasal bone, free beta-hCG and PAPP-A at 11 to 14 weeks. *Prenat Diagn.* 2003 Apr;23(4):306-10.
10. Maiz N, Valencia C, Kagan KO, Wright D, Nicolaides KH. Ductus venosus Doppler in screening for trisomies 21, 18 and 13 and Turner syndrome at 11-13 weeks of gestation. *Ultrasound Obstet Gynecol.* 2009 May;33(5):512-7.
11. Kagan KO, Valencia C, Livanos P, Wright D, Nicolaides KH. Tricuspid regurgitation in screening for trisomies 21, 18 and 13 and Turner syndrome at 11+0 to 13+6 weeks of gestation. *Ultrasound Obstet Gynecol.* 2009 Jan;33(1):18-22.
12. Sonek J, Borenstein M, Dagklis T, Persico N, Nicolaides KH. Frontomaxillary facial angle in fetuses with trisomy 21 at 11-13(6) weeks. *Am J Obstet Gynecol.* 2007 Mar;196(3):271.e1-4.
13. Chaoui R, Benoit B, Mitkowska-Wozniak H, Heling KS, Nicolaides KH. Assessment of intracranial translucency (IT) in the detection of spina bifida at the 11-13-week scan. *Ultrasound Obstet Gynecol* 2009; 34: 249-252.
14. Ozkaya O, Sezik M, Ozbasar D, Kaya H. Abnormal ductus venosus flow and tricuspid regurgitation at 11-14 weeks' gestation have high positive predictive values for increased risk in first-trimester combined screening test: results of a pilot study. *Taiwan J Obstet Gynecol.* 2010 Jun;49(2):145-50.
15. Bethune M. Literature review and suggested protocol for managing ultrasound soft markers for Down syndrome: thickened nuchal fold, echogenic bowel, shortened femur, shortened humerus, pyelectasis and absent or hypoplastic nasal bone. *Australas Radiol.* 2007 Jun;51(3):218-25
16. Nicolaides KH. Screening for fetal aneuploidies at 11 to 13 weeks. *Prenat Diagn.* 2011 Jan;31(1):7-15. doi: 10.1002/pd.2637.
17. Valinen Y, Marttala J, Sahravand M, Jarvela I, Honkasalo T, Laitinen P, Ryyänen M. Adding ADAM12 in risk calculation program does not improve the detection rate of trisomies 18 and 13 in first trimester screening. *J Matern Fetal Neonatal Med.* 2012 Feb 6
18. Valinen Y, Laitinen P, Ranta J, Ignatius J, Jarvela I, Ryyänen M. Effect of a new marker, ADAM12, on Down risk figures in first trimester screening. *J Matern Fetal Neonatal Med.* 2009 Jul;22(7):602-7.

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Bibliografia

1. Snijders RJ, Noble P, Sebire N, Souka A, Nicolaides KH. UK multicentre project on assessment of risk of trisomy 21 by maternal age and fetal nuchal-translucency thickness at 10-14 weeks of gestation. Fetal Medicine Foundation First Trimester Screening Group. *Lancet*. 1998 Aug 1;352(9125):343-6.
2. Jou HJ, Wu SC, Li TC, Hsu HC, Tzeng CY, Hsieh FJ. Relationship between fetal nuchal translucency and crown-rump length in an Asian population. *Ultrasound Obstet Gynecol*. 2001 Feb;17(2):111-4.
3. Salomon LJ, Chalouhi GE, Bernard JP, Ville Y; Société française pour l'amélioration des pratiques échographiques (SFAPE). [Nuchal translucency thickness at 11-14 weeks of gestation: French charts and equations. *J Gynecol Obstet Biol Reprod (Paris)*. 2009 Dec;38(8):635-41.
4. Kor-Anantakul O, Suntharasaj T, Suwanrath C, Chanprapaph P, Sirichotiyakul S, Ratanasiri T, Vuthiwong C, Geater A. Distribution of normal nuchal translucency thickness: a multicenter study in Thailand. *Gynecol Obstet Invest*. 2011;71(2):124-8.
5. Markov D, Chernev T, Dimitrova V, Mazneŭkova V. First trimester nuchal translucency thickness in normal fetuses in the Bulgarian population. *Akush Ginekol (Sofia)*. 2005;44(5):3-8.
6. Chung JH, Yang JH, Song MJ, Cho JY, Lee YH, Park SY, Moon MJ, Lim HJ, Choi JS, Kim JO, Shin JS, Ahn HK, Han JY, Kim MY, Choi KH, Ryu HM. The distribution of fetal nuchal translucency thickness in normal Korean fetuses. *J Korean Med Sci*. 2004 Feb;19(1):32-6.
7. Zoppi MA, Ibba RM, Floris M, Monni G. Fetal nuchal translucency screening in 12495 pregnancies in Sardinia. *Ultrasound Obstet Gynecol*. 2001 Dec;18(6):649-51.
8. Spencer K, Bindra R, Nix AB, Heath V, Nicolaides KH. Delta-NT or NT MoM: which is the most appropriate method for calculating accurate patient-specific risks for trisomy 21 in the first trimester? *Ultrasound Obstet Gynecol*. 2003 Aug;22(2):142-8.
9. Wright D, Kagan KO, Molina FS, Gazzoni A, Nicolaides KH. A mixture model of nuchal translucency thickness in screening for chromosomal defects. *Ultrasound Obstet Gynecol*. 2008 Apr;31(4):376-83.
10. Bestwick JP, Huttly WJ, Wald NJ. Evaluation of a proposed mixture model to specify the distributions of nuchal translucency measurements in antenatal screening for Down's syndrome. *J Med Screen*. 2010;17(1):13-8.
11. Novelli A, Grati FR, Ballarati L, Bernardini L, Bizzoco D, Camurri L, Casalone R, Cardarelli L, Cavalli P, Ciccone R, Clementi M, Dalprà L, Gentile M, Gelli G, Grammatico P, Malacarne M, Nardone AM, Pecile V, Simoni G, Zuffardi O, Giardino D. Microarray application in prenatal diagnosis: a position statement from the cytogenetics working group of the Italian Society of Human Genetics (SIGU), November 2011. *Ultrasound Obstet Gynecol*. 2012 Apr;39(4):384-8.

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Bibliografia

1. Filly RA, Benacerraf BR, Nyberg DA, Hobbins JC. Choroid plexus cyst and echogenic intracardiac focus in women at low risk for chromosomal anomalies. *J Ultrasound Med.* 2004 Apr;23(4):447-9.
2. Bethune M. Literature review and suggested protocol for managing ultrasound soft markers for Down syndrome: thickened nuchal fold, echogenic bowel, shortened femur, shortened humerus, pyelectasis and absent or hypoplastic nasal bone. *Australas Radiol.* 2007 Jun;51(3):218-25

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Figura 3 Plica nucale aumentata di 6.2 mm (Bcolor)

Bibliografia

1. Benacerraf BR, Neuberg D, Bromley B, Frigoletto FD Jr. Sonographic scoring index for prenatal detection of chromosomal abnormalities. *J Ultrasound Med.* 1992 Sep;11(9):449-58.
2. Nyberg DA, Souter VL, El-Bastawissi A, Young S, Luthhardt F, Luthy DA. Isolated sonographic markers for detection of fetal Down syndrome in the second trimester of pregnancy. *J Ultrasound Med.* 2001 Oct;20(10):1053-63.
3. Raniga S, Desai PD, Parikh H. Ultrasonographic soft markers of aneuploidy in second trimester: are we lost? *MedGenMed.* 2006 Jan 11;8(1):9.
4. Benacerraf BR. The role of the second trimester genetic sonogram in screening for fetal Down syndrome. *Semin Perinatol.* 2005 Dec;29(6):386-94.
5. Van den Hof MC, Wilson RD; Diagnostic Imaging Committee, Society of Obstetricians and Gynaecologists of Canada; Genetics Committee, Society of Obstetricians and Gynaecologists of Canada. Fetal soft markers in obstetric ultrasound. *J Obstet Gynaecol Can.* 2005 Jun;27(6):592-636.
6. Cho JY, Kim KW, Lee YH, Toi A. Measurement of nuchal skin fold thickness in the second trimester: influence of imaging angle and fetal presentation. *Ultrasound Obstet Gynecol.* 2005 Mar;25(3):253-7.
7. Smith-Bindman R, Hosmer W, Feldstein VA, Deeks JJ, Goldberg JD. Second-trimester ultrasound to detect fetuses with Down syndrome: a meta-analysis. *JAMA.* 2001 Feb 28;285(8):1044-55.

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Figura 1 Focus iperecogeno cardiaco

Bibliografia

1. Anderson N, Jyoti R. Relationship of isolated fetal intracardiac echogenic focus to trisomy 21 at the mid-trimester sonogram in women younger than 35 years. *Ultrasound Obstet Gynecol.* 2003 Apr;21(4):354-8.
2. De Vore G.R. et al; OP01.16, 16th World Congress on Ultrasonod in Obstet and Gynecol London 2006
3. Smith-Bindman R, Hosmer W, Feldstein VA, Deeks JJ, Goldberg JD. Second-trimester ultrasound to detect fetuses with Down syndrome: a meta-analysis. *JAMA.* 2001 Feb 28;285(8):1044-55.
4. Coco C, Jeanty P, Jeanty C. An isolated echogenic heart focus is not an indication for amniocentesis in 12,672 unselected patients. *J Ultrasound Med.* 2004 Apr;23(4):489-96.

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Bibliografia

1. Samánek M, Vorísková M. Congenital heart disease among 815,569 children born between 1980 and 1990 and their 15-year survival: a prospective Bohemia survival study. *Pediatr Cardiol.* 1999 Nov-Dec;20(6):411-7.
2. Meberg A, Otterstad JE, Frøland G, Sørland S, Nitter-Hauge S.. Increasing incidence of ventricular septal defects caused by improved detection rate. *Acta Paediatr.* 1994 Jun;83(6):653-7.
3. Paladini D, Palmieri S, Lamberti A, Teodoro A, Martinelli P, Nappi C. Characterization and natural history of ventricular septal defects in the fetus. *Ultrasound Obstet Gynecol.* 2000 Aug;16(2):118-22.
4. Axt-Flidner R, Schwarze A, Smrcek J, Germer U, Krapp M, Gembruch U. Isolated ventricular septal defects detected by color Doppler imaging: evolution during fetal and first year of postnatal life. *Ultrasound Obstet Gynecol.* 2006 Mar;27(3):266-73.
5. Glen S, Burns J, Bloomfield P. Prevalence and development of additional cardiac abnormalities in 1448 patients with congenital ventricular septal defects. *Heart.* 2004 Nov;90(11):1321-5.
6. Li H, Wei J, Ma Y, Shang T. Prenatal diagnosis of congenital fetal heart abnormalities and clinical analysis. *J Zhejiang Univ Sci B.* 2005 Sep;6(9):903-6.
7. Tennstedt C, Chaoui R, Körner H, Dietel M. Spectrum of congenital heart defects and extracardiac malformations associated with chromosomal abnormalities: results of a seven year necropsy study. *Heart.* 1999 Jul;82(1):34-9.
8. Hyett J, Moscoso G, Nicolaides K. Abnormalities of the heart and great arteries in first trimester chromosomally abnormal fetuses. *Am J Med Genet.* 1997 Mar 17;69(2):207-16.
9. Hyett J, Perdu M, Sharland G, Sniijders R, Nicolaides KH. Using fetal nuchal translucency to screen for major congenital cardiac defects at 10-14 weeks of gestation: population based cohort study. *BMJ.* 1999 Jan 9;318(7176):81-5.
10. Moyano D, Huggon IC, Allan LD. Fetal echocardiography in trisomy 18. *Arch Dis Child Fetal Neonatal Ed.* 2005 Nov;90(6):F520-2.
11. McElhinney DB, Driscoll DA, Levin ER, Jawad AF, Emanuel BS, Goldmuntz E. Chromosome 22q11 deletion in patients with ventricular septal defect: frequency and associated cardiovascular anomalies. *Pediatrics.* 2003 Dec;112(6 Pt 1):e472.
12. Tennstedt C, Chaoui R, Körner H, Dietel M. Spectrum of congenital heart defects and extracardiac malformations associated with chromosomal abnormalities: results of a seven year necropsy study. *Heart.* 1999 Jul;82(1):34-9.
13. Hyett J, Moscoso G, Nicolaides K. Abnormalities of the heart and great arteries in first trimester chromosomally abnormal fetuses. *Am J Med Genet.* 1997 Mar 17;69(2):207-16.
14. Tennstedt C, Chaoui R, Körner H, Dietel M. Spectrum of congenital heart defects and extracardiac malformations associated with chromosomal abnormalities: results of a seven year necropsy study. *Heart.* 1999 Jul;82(1):34-9.
15. Paladini D, Tartaglione A, Agangi A, Teodoro A, Forleo F, Borghese A, Martinelli P. The association between congenital heart disease and Down syndrome in prenatal life. *Ultrasound Obstet Gynecol.* 2000 Feb;15(2):104-8.
16. DeVore GR. Trisomy 21: 91% detection rate using second-trimester ultrasound markers. *Ultrasound Obstet Gynecol.* 2000 Aug;16(2):133-41.

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Bibliografia

1. Nyberg DA, Resta RG, Luthy DA, Hickok DE, Mahony BS, Hirsch JH. Prenatal sonographic findings of Down syndrome: review of 94 cases. *Obstet Gynecol.* 1990 Sep;76(3 Pt 1):370-7.
2. Benacerraf BR, Pober BR, Sanders SP. Accuracy of fetal echocardiography. *Radiology.* 1987 Dec;165(3):847-9.
3. DeVore GR. Trisomy 21: 91% detection rate using second-trimester ultrasound markers. *Ultrasound Obstet Gynecol.* 2000 Aug;16(2):133-41.
4. Huisman TW, Stewart PA, Stijnen T, Wladimiroff JW. Doppler flow velocity waveforms in late first- and early second-trimester fetuses: reproducibility of waveform recordings. *Ultrasound Obstet Gynecol.* 1993 Jul 1;3(4):260-3.
5. Messing B, Porat S, Imbar T, Valsky DV, Anteby EY, Yagel S. Mild tricuspid regurgitation: a benign fetal finding at various stages of pregnancy. *Ultrasound Obstet Gynecol.* 2005 Nov;26(6):606-9; discussion 610.
6. Huggon IC, DeFigueiredo DB, Allan LD. Tricuspid regurgitation in the diagnosis of chromosomal anomalies in the fetus at 11-14 weeks of gestation. *Heart.* 2003 Sep;89(9):1071-3.
7. Kagan KO, Valencia C, Livanos P, Wright D, Nicolaides KH. Tricuspid regurgitation in screening for trisomies 21, 18 and 13 and Turner syndrome at 11+0 to 13+6 weeks of gestation. *Ultrasound Obstet Gynecol.* 2009 Jan;33(1):18-22.

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
Figura 1 Osso nasale

Bibliografia

1. Cusick W, Provenzano J, Sullivan CA, Gallousis FM, Rodis JF. Fetal nasal bone length in euploid and aneuploid fetuses between 11 and 20 weeks' gestation: a prospective study. *J Ultrasound Med* 2004; 23: 1327-1333.
2. Bromley B, Lieberman E, Shipp TD, Benacerraf BR. Fetal nasal bone length: a marker for Down syndrome in the second trimester. *J Ultrasound Med* 2002; 21: 1387-1394.
3. Viora E, Errante G, Sciarrone A, Bastonero S, Masturzo B, Martiny G, Campogrande M. Fetal nasal bone and trisomy 21 in the second trimester. *Prenat Diagn* 2005; 25: 511-515.
4. Odibo AO, Sehdev HM, Dunn L, McDonald R, Macones GA. The association between fetal nasal bone hypoplasia and aneuploidy. *Obstet Gynecol* 2004; 104: 1229-1233.
5. Cicero S, Sonek JD, McKenna DS, Croom CS, Johnson L, Nicolaides KH. Nasal bone hypoplasia in trisomy 21 at 15-22 weeks' gestation. *Ultrasound Obstet Gynecol* 2003; 21:15-18.
6. Sonek JD, Nicolaides KH. Prenatal ultrasonographic diagnosis of nasal bone abnormalities in three fetuses with Down syndrome. *Am J Obstet Gynecol* 2002; 186: 139-141.

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Bibliografia

1. Ouahba J, Luton D, Vuillard E, Garel C, Gressens P, Blanc N, Elmaleh M, Evrard P, Oury JF. Prenatal isolated mild ventriculomegaly: outcome in 167 cases. *BJOG*. 2006 Sep;113(9):1072-9.
2. Papp C, Beke A, Ban Z, Szigeti Z, Toth-Pal E, Papp Z. Prenatal diagnosis of trisomy 13: analysis of 28 cases. *J Ultrasound Med*. 2006 Apr;25(4):429-35.
3. Nyberg DA, Souter VL. Sonographic markers of fetal trisomies: second trimester. *J Ultrasound Med*. 2001 Jun;20(6):655-74.
4. Bronsteen R, Lee W, Vettrano IM, Huang R, Comstock CH. Second-trimester sonography and trisomy 18: the significance of isolated choroid plexus cysts after an examination that includes the fetal hands. *J Ultrasound Med*. 2004 Feb;23(2):241-5.
5. Ghidini A, Strobelt N, Locatelli A, Mariani E, Piccoli MG, Vergani P. Isolated fetal choroid plexus cysts: role of ultrasonography in establishment of the risk of trisomy 18. *Am J Obstet Gynecol*. 2000 Apr;182(4):972-7.
6. Whitlow BJ, Lazanakis ML, Kadir RA, Chatzipapas I, Economides DL. The significance of choroid plexus cysts, echogenic heart foci and renal pyelectasis in the first trimester. *Ultrasound Obstet Gynecol*. 1998 Dec;12(6):385-90.
7. Coco C, Jeanty P. Isolated fetal pyelectasis and chromosomal abnormalities. *Am J Obstet Gynecol*. 2005 Sep;193(3 Pt 1):732-8.
8. Ecker JL, Shipp TD, Bromley B, Benacerraf B. The sonographic diagnosis of Dandy-Walker and Dandy-Walker variant: associated findings and outcomes. *Prenat Diagn*. 2000 Apr;20(4):328-32.
9. Blaicher W, Ulm B, Ulm MR, Hengstschläger M, Deutinger J, Bernaschek G. Dandy-Walker malformation as sonographic marker for fetal triploidy. *Ultraschall Med*. 2002 Apr;23(2):129-33.
10. Chen CP. Aplasia and duplication of the thumb and facial clefts associated with fetal trisomy 18. *Am J Med Genet A*. 2006 May 1;140(9):960-3.
11. Perrotin F, de Poncheville LM, Marret H, Paillet C, Lansac J, Body G. Chromosomal defects and associated malformations in fetal cleft lip with or without cleft palate. *Eur J Obstet Gynecol Reprod Biol*. 2001 Nov;99(1):19-24.
12. Bollmann R, Kalache K, Mau H, Chaoui R, Tennstedt C. Associated malformations and chromosomal defects in congenital diaphragmatic hernia. *Fetal Diagn Ther*. 1995 Jan-Feb;10(1):52-9.
13. Groves R, Sunderajan L, Khan AR, Parikh D, Brain J, Samuel M. Congenital anomalies are commonly associated with exomphalos minor. *J Pediatr Surg*. 2006 Feb;41(2):358-61.
14. Oliveira EA, Cabral AC, Pereira AK, Machado IN, Diniz JS, Lana AM, Cunha LK. Outcome of fetal urinary tract anomalies associated with multiple malformations and chromosomal abnormalities. *Prenat Diagn*. 2001 Feb;21(2):129-34.
15. Favre R, Köhler M, Gasser B, Müller F, Nisand I. Early fetal megacystis between 11 and 15 weeks of gestation. *Ultrasound Obstet Gynecol*. 1999 Dec;14(6):402-6.
16. Chen CP, Chern SR, Lee CC, Town DD, Chen WL, Wang W. Bilateral renal agenesis and fetal ascites in association with partial trisomy 13 and partial trisomy 16 due to a 3:1 segregation of maternal reciprocal translocation t(13;16)(q12.3; p13.2). *Prenat Diagn*. 1999 Aug;19(8):783-6.
17. Cocchi G, Magnani C, Morini MS, Garani GP, Milan M, Calzolari E. Urinary tract abnormalities (UTA) and associated malformations: data of the Emilia-Romagna Registry. IMER Group. Emilia-Romagna Registry on Congenital Malformations. *Eur J Epidemiol*. 1996 Oct;12(5):493-7.
18. Bornstein E, Barnhard Y, Donnenfeld A, Ferber A, Divon MY. Fetal pyelectasis: does fetal gender modify the risk of major trisomies? *Obstet Gynecol*. 2006 Apr;107(4):877-9.
19. Wiesel A, Queisser-Luft A, Clementi M, Bianca S, Stoll C; EUROSCAN Study Group. Prenatal detection of congenital renal malformations by fetal ultrasonographic examination: an analysis of 709,030 births in 12 European countries. *Eur J Med Genet*. 2005 Apr-Jun;48(2):131-44.
20. Smilari P, La Spina M, Salvo V, Romeo MG, Sanges G. Esophageal atresia and malformative association. Clinical contribution. *Minerva Pediatr*. 2005 Oct;57(5):289-96.
21. Heinen FL, Vallone P, Elmo G. Esophageal diverticulum in an infant with Down's syndrome and type III esophageal atresia. *J Pediatr Surg*. 2003 Apr;38(4):E9.
22. Martínez-Frías ML, Bermejo E. Major congenital malformations in Down syndrome. *Am J Med Genet*. 1997 Nov 28;73(1):91-2.
23. Pameijer CR, Hubbard AM, Coleman B, Flake AW. Combined pure esophageal atresia, duodenal atresia, biliary atresia, and pancreatic ductal atresia: prenatal diagnostic features and review of the literature. *J Pediatr Surg*. 2000 May;35(5):745-7.

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Modificato da: Wald NJ et al. First and second trimester antenatal screening or Down's syndrome: the results of the Serum, Urine and Ultrasound Screening Study (SURUSS).

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Bibliografia

1. Wright D, Bradbury I, Benn P, Cuckle H, Ritchie K. Contingent screening for Down syndrome is an efficient alternative to non-disclosure sequential screening. *Prenat Diagn.* 2004 Oct;24(10):762-6.
2. Filkins K, Koos BJ. Ultrasound and fetal diagnosis. *Curr Opin Obstet Gynecol.* 2005 Apr;17(2):185-95.
3. Cleary-Goldman J, Morgan MA, Malone FD, Robinson JN, D'Alton ME, Schulkin J. Screening for Down syndrome: practice patterns and knowledge of obstetricians and gynecologists. *Obstet Gynecol.* 2006 Jan;107(1):11-7.